

1. Sterilization system

Designation	Autoclave
Visualization	
Manufacture	Belimed Sauter AG, CH-8583 Sulgen
Factory number	6-6-12 VS2
Internal name	Sterylizator 2
Year of build	2007
Containing	1 trolley / 2 layers
operating steam pressure	2,7bar
Kind of process	Steam

The periferie is described in the TD-Study. See report 2404950

1.1 Perifery

Pumps / vaveles / sensors



Figure 1: Steam generator



Figure 2: compressed air, water



Figure 3: Circulation pump



Figure 4: Water valve



Figure 5: Temperatur sensor

All parts are located on the right-hand side as seen from the front.




Figure 15: Visualization of cold spot (blue)